\equiv



MedAssist Live - Caretaker

MedAssist Live connects a caretaker to a coach, anytime, anywhere.







A Coach instructs a Caretaker

A Caretaker works with a patient/client

1. Log in

Step	Action	
1	Check that you have a 4G or WIFI connection	(îŗ
2	Download from your App Store the MedAssist Live app Open the Live app	Live

2. Connect screen

2.1 Collaboration without case management

In this option the caretaker clicks on a link, that was sent by a coach.

Step	Action		
1	Click on the link that you received via email, WhatsApp or SMS. The message should look like this:	Please join me in MedAssist Live. https://medassistlive.page.link/kfSU	
2	MedAssist Live starts up and you enter the call immediately		
	 You open the app WITHOUT logging in. No user ID required. Wait for a Two Way Augmented Reality (TWAR) session invitation to appear on the screen. Click on Menu options are limited when connecting through a link. 	You are invited to join TWAR session #278839 Cancel GO	

2.2 Collaboration with case management

The purpose of the connect screen is to enable the caretaker choose a coach to connect with and to call the coach, when available.

! Please view the user manual for instructions on how to use case management.

Step	Action				
1	Enter User ID and password				
	! Only the first tim before moving in application.	ne, agree to the Terms of Use to the Main Menu of the			
2	Case Management enabled If you have case management enabled the following menu option appear:				
	Caretaker Option				
	Connect	Connect with a coach using this menu option.			
	Cases	A Case can be created and maintained with information about the patient, the incident, symptoms, advice given.			
	Patients A Patient record can be created, changed and deleted.				
	Profile Screen	(Right upper side of the screen): User details,			

Access to documentation, team information.

version information, ability to log out.



	The screen looks like this:	Caretaker	3. Co This scr	Dllaboration screen reen connects the caretaker and a coach to work in the same virtual reality.
		Connect Do you need medical assistance? Connect Cases © View all cases © Add new case	Step 1	 Action The caretaker should direct the device with the front camera towards the patient or client. The screen should be in Landscape
		Patients [©] View all patients G Add new patient C C C C C C C C C C	2	mode, otherwise a warning message will appear. When the connection has been successfully established, the caretaker and coach can
	Case Management disabled If you do not have case management you will see the connect screen.	Caretaker Careta		 Collaborate Review Case information Flip the camera to have face-to-face communication Switch their microphone on/off Exchange chat messages End the call
			4. Re You car <u>https://</u>	eferences n find the user manual and other documentation via this link : /medassist.online/liveapp/docs